

## Pulmo-Aide® Compressor Nebulizer System

### The Gold Standard for High Performance Medication Delivery

The DeVilbiss Pulmo-Aide has been providing **effective aerosol treatments for more than thirty years**. Modern styling, durability and ease of use make this unit **perfect for in-home asthma treatments**. The Pulmo-Aide comes ready to use with all the attachments - filter, tubing, and nebulizer.

The **Pulmo-Aide is the only approved compressor for use with TOBI** (tobramycin solution for inhalation). The DeVilbiss Pulmo-Aide is unmatched in quality, value and performance.



Pulmo-Aide Compressor  
Nebulizer System  
Model 5650D

### Performance Features

- Maximum compressor pressure – 35 psig or greater
- Compressor free air flow – 9 lpm or greater
- 5-year limited warranty
- Longest lasting compressor on the market with an average product life of 10 years
- Unmatched quality and durability - Great for nursing homes and long-term users

### Design Features

- Metal frame motor construction
- Unique 4-pole motor design extends life of the unit
- Large built-in compartment for storing nebulizers and tubing
- Built-in nebulizer holder is a convenient place to store your nebulizer if your treatment is interrupted.

## Available Models

<b>5650D</b>	Pulmo-Aide Compressor Nebulizer System 115VAC, 2 prong plug
<b>5650H</b>	Pulmo-Aide Compressor Nebulizer System 230VAC, 3 prong plug
<b>5650C</b>	Pulmo-Aide Compressor Nebulizer System CSA approved

## Parts and Accessories

<b>DEVILBISS PART #</b>	<b>PRODUCT DESCRIPTION</b>
<b>800D</b>	Reusable Nebulizer – blue (1 each)
<b>800C</b>	Reusable Nebulizer – purple/green (1 each)
<b>4650D-621</b>	Disposable Nebulizer (50 pack)
<b>AAM12</b>	Adult Aerosol Mask (12 pack)
<b>AAM20</b>	Adult Aerosol Mask (50 pack)
<b>PAM10</b>	Pediatric Aerosol Mask (50 pack)
<b>DLI050</b>	Pediatric Children's Character Mask (50 pack)
<b>5650D-601</b>	Compressor Filter (package of 25)
<b>5650D-602</b>	Compressor Filter (package of 5)
<b>5650D-612</b>	Filter Kits (6 each – caps and filters)
<b>5650D-611</b>	Gray Hose Tubing Connector (package of 6)

## 800 Reusable Jet Nebulizer

- MMAD 3.6 – 4 microns
- Approximately 75% of particles are below 5 microns in size
- Reduces treatment time nearly in half – 3cc treatment in 7 – 9 minutes
- 6cc capacity
- Dishwasher safe to promote frequent cleaning
- Ergonomic design fits comfortably in the hand – even with limited functionality
- 4 components for ease of assembly/disassembly
- Comfortable mouthpiece design
- Standard 22mm fitting, no adapters required
- Can be used inline with a ventilator
- 1 year warranty
- Tested for use with Albuterol, Pulmicort (budesonidel, Ipratropium bromide)
- Suitable for use with most common liquid medications

### Technical Specifications

<b>Built to meet the requirements of HCPCS code:</b> E0570	
<b>Product Dimensions:</b>	10.5H X 6.3W X 10.1D
<b>Product Weight:</b>	7.1 lbs.
<b>Maximum Air Pressure:</b>	30+ psig
<b>Maximum Liter Flow:</b>	9+ lpm
<b>MMAD:</b>	5 microns or less
<b>Sound:</b>	51 dBA
<b>Filter Replacement:</b>	No tools required
<b>Warranty:</b>	5 Years
<b>Safety Approvals:</b>	UL1431
<b>Motor/Pump Assembly:</b>	Diaphragm, Primarily Metal Parts
<b>Neb Storage:</b>	Yes
<b>Master Carton Quantity:</b>	4



### The HomeCare Medical Ltd.

401-403 Wing Ming Industrial Centre,  
15 Cheung Yue Street,  
Lai Chi Kok, Kowloon  
Phone: 2402-2188 Fax: 2411-2323