



spiro lab[®] III

diagnostic spirometer
with oximetry option

New!
3rd generation

Fast
Simple
Durable
Accurate
Powerful
Innovative

- Spirometer with 6,000 test memory
- Oximeter with 1,000 hours of recording
- Available with disposable or reusable turbine flowmeter



Remote diagnosis via web also with Bluetooth[®] connection



Telemedicine
option

SpO2
option



PC software included



Quality Instruments

FDA	ATS	ISO	ISO
Approved	Standard	9001-2000	13485



Quality Instruments

spirolab[®] III
complies with
the new ATS
standard


Spirometer

- Records the best 3 trials
- Up to 8 tests on one screen
- Internal temperature sensor for BTPS conversion



SpO2 Option

- Short or long term (overnight) SpO2 and Pulse Rate measurements
- ODI, NOD, T90%, T89%, T88%, T87%
- Sleep oximetry with desaturation events
- Oximetry during exercise test



Exclusive pediatric incentive system

spirolab[®] III

diagnostic spirometer with oximetry option

Standard device includes:
spirolab[®] III
winspiroPRO CD
Carrying case

FVC, VC with breathing pattern plus MVV tests with real time curves

- High resolution color screen
- Fast and silent thermal printer
- Digital turbine flowmeter with guaranteed accuracy in all environmental condition
- PRE-POST bronchodilator comparison
- Selectable language and predicted values
- Connectivity: USB, **Bluetooth[®]** and RS232



Two different flowmeters available

FlowMir disposable turbine for single patient use

- High accuracy
- Easy to replace
- Cost Effective
- Individually tested and calibrated
- Available in box of 100 pieces
- No maintenance – No filter – No problem
- Hygiene 100% guaranteed by single packaging



MIR Exclusive Patent

Reusable turbine for long term operation

- High accuracy
- Long term stability
- Easy to clean or sanitized



MIR digital reusable and disposable turbines are developed in full compliance with ATS standards and guarantee accuracy in all environmental conditions



