

Choose your **GO₂ Achieve** Accessories



GO2CC

Black neoprene carrying case with belt loop



GO2L

Black 20" Lanyard



GO2R

Retractable Clip Holder



Three colors to choose from



GO₂ Achieve
Fingertip Pulse Oximeter

Specifications

Oxygen Saturation Display Range	0% to 100
Pulse Rate Display Range	18 to 321 pulses per minute
Oxygen Saturation Declared Accuracy Range (Arms*)	70-100% SpO ₂ ± 2 digits
Low Perfusion Oxygen Saturation Declared Accuracy Range (Arms*)	70-100% SpO ₂ ± 2 digits
Pulse Rate Declared Accuracy Range (Arms*)	20-250 BPM ± 3 digits
Low Perfusion Pulse Rate Declared Accuracy Range (Arms*)	40-240 ± 3 digits

*± 1 Arms represents approximately 68% of measurements.

Temperature (Operating)	+41° to +104°F (5°C to +40°C)
Storage/Transportation	-22° to +158°F (-30° to +70°C)
Device temperature will not exceed 41°C as measured in a controlled environment.	
Humidity (Operating)	10% to 90% relative humidity, non-condensing
Storage/Transportation	10% to 95% relative humidity, non-condensing
Operating Altitude	Up to 40,000 feet / 12,192 meters
Battery Life (Continuous):	Approximately 2,400 spot checks based on ~21 hours of operation using one AAA-size alkaline battery, calculated at 30 seconds per use.
Battery Life (Storage)	6 months minimum

2 year warranty.

Backlit LCD for low light viewing.

Specifications are subject to change without notice.



Nonin Medical, Inc.

13700 1st Avenue North • Plymouth MN • 55441-5443 • USA
Phone: +1.763.553.9968 • Toll Free: 877.577.2635 • Fax: +1.763.577.5521
E-mail: info@nonin.com



www.go2nonin.com



Authorized Distributor:

©2009 Nonin Medical, Inc.
All trademarks are the property of Nonin Medical, Inc. unless otherwise noted.
7070-000-03



Go Further with **GO₂ Achieve**

Reliable oxygen saturation and pulse rate monitoring for any environment




Fingertip Pulse Oximeter

GO₂ Achieve

Nonin Medical, the inventor of fingertip pulse oximetry, introduces the **GO₂ Achieve** fingertip pulse oximeter. The affordable **GO₂ Achieve** provides you with the knowledge of reliable oxygen saturation and pulse rate values without compromising the quality and accuracy that have made Nonin the most trusted name in fingertip pulse oximetry.

Live life on the GO!



GO₂ Achieve



Get Fit with GO₂ Achieve

The highly durable **GO₂ Achieve** fingertip oximeter is an ideal oxygen and pulse rate monitoring companion while you exercise. Compact and easy-to-use, simply place **GO₂ Achieve** on your finger for instant accurate oxygen saturation and pulse rate feedback.



Fly with Peace of Mind

When flying at high altitudes, hypoxia can impair cognitive function beginning at 5,000 feet. For instant and reliable results, the **GO₂ Achieve** easily provides the necessary information for passengers and pilots to check for the onset of in-flight hypoxia.



Blood Oxygen Saturation By Altitude

Altitude (feet)	Arterial O ₂ Saturation Without Supplemental O ₂	Atmospheric Pressure (mmHg)
0	96%	760
5,000	95%	632
7,500	93%	575
10,000	89%	523
12,500	87%	474
14,000	83%	446
16,500	77%	403
20,000	65%	349
25,000	Below 60%	282

"The Pilot: An Air Breathing Mammal," Mehler, Stanley R. MD, Human Factors Bulletin, Flight Safety Foundation, 1981.



Climb Higher

The accurate and durable **GO₂ Achieve** is the most dependable way to monitor for Acute Mountain Sickness. Nonin's products are recognized by the American Alpine Institute and have been to the summit of many of the world's highest peaks including Denali (20,320') and Everest (29,035').

Quality and Support DO Matter

With so many throwaway products available these days, it's important to choose a product that not only is consistently accurate, but is also built to last. At Nonin, we are dedicated to providing the highest quality products available at a great value.

Superior Service and Support

Nonin will be here long after the sale and our knowledgeable customer service representatives are always available to answer your questions.

A Message for All

Not All Oximeters are Created Equal!

- Nonin welcomes scrutiny of our independent testing results and they are available for all to see! Check out why Nonin is the most trusted name in fingertip oximetry at www.go2nonin.com.
- Nonin's quality standards require every **GO₂ Achieve** be able to withstand 25 drops from one meter onto a hard surface. Try it with a **GO₂ Achieve** imitator and be prepared to pick up the pieces!
- Life can get messy — **GO₂ Achieve** exceeds IP33 water ingress testing.
- The **GO₂ Achieve** is NOT a throw-away product. It is built to last and covered by our 2-year warranty. Nonin is a company that will be here tomorrow.
- At Nonin, you'll speak with a real person right away. The voice you hear represents the pride in our products. We go above and beyond to support our satisfied customers.