

DeVilbiss Traveler® Portable Compressor Nebulizer System

The **perfect** travel companion



DeVilbiss Traveler Portable Compressor Nebulizer Unit with battery attached Model 6910P-DR

The **DeVilbiss Traveler Portable Compressor Nebulizer System** is designed to help active customers get more out of their aerosol therapy, and **more out of life**. It's never been easier or more convenient to nebulize outside of the home.

The DeVilbiss Traveler includes compressor unit, **disposable** and reusable nebulizers, **five extra filters**, tubing adapter, adult and **pediatric** masks, AC adapter, DC adapter, carrying case, instruction guide and **instructional DVD**, as well as a **full-year** battery warranty. Units can be ordered with **optional battery**.

Performance Features

- Delivers a **consistent and concentrated treatment** of small particles that penetrate deeply into the lungs, delivering medication where it's needed most
- **Fast treatment time** (approximately 10 minutes) with DeVilbiss 900 Series Reusable Nebulizer encourages patient compliance
- Can be used with **any manufacturer's nebulizer** – treatment times will vary
- **45-minute battery life** far outlasts treatment time

Design Features

- Easy, one-touch operation
- Compact and lightweight
- AC adapter and DC car adapter
- One-year battery and three-year compressor warranty

Model 6910D-DR Includes:

- Compressor unit
- Reusable nebulizer with tubing
- Disposable nebulizer
- 5 replacement air filters and filter door
- Tubing adapter
- Adult and pediatric masks
- AC adapter
- DC car adapter
- Easy tote carrying case
- Instructional manual and DVD

Model 6910P-DR Includes:

- All of the above items, plus a rechargeable battery pack

Free Consumer Marketing Program to Generate Sales

- Comprehensive Traveler marketing kit includes counter card with reminder tear pad, product flyer (staple to prescription bags) and patient postcard mailers to generate business



DeVilbiss Traveler Portable Compressor Nebulizer System Model 6910P-DR

Your nebulizer may be different than shown.

Accessories

6910D-601	Rechargeable Battery
6910D-602	AC Adapter
6910D-603	DC Car Adapter
6910D-605	Air-Inlet Filter (5), Filter Cover (1), & Neb Bracket (1)
6910D-607	Carrying Case
900D-610	Reusable Nebulizer (10/case)
3655D-621	Disposable Nebulizer (50/case)
932168	Instructional DVD
932167	Instructional VHS Video (NTSC format)
AAM12	Adult Aerosol Mask (12/case)
PAM10	Pediatric Aerosol Mask (50/case)
DL1050	Pediatric Character Mask (50/case)



Technical Specifications

MMAD	5 microns	Sound Level	54 dBA
Size w/o Battery	3.34" W x 2.2" H x 3.8" D	Max. Pressure	30 psig or greater
Size w/ Battery	3.34" W x 2.2" H x 5.3" D	Free Air Flow	6 lpm or greater
Weight w/o Battery	.8 lbs	Intermittent Use	30 minutes on/30 minutes off
Weight w/ Battery	1.5 lbs	Master Carton Size	8 units
Electrical	100-240V, 50-60Hz, 12VDC, 1.5A	Battery Life	45 minutes

The HomeCare Medical Ltd.

401-403 Wing Ming Industrial Centre,
15 Cheung Yue Street,
Lai Chi Kok, Kowloon
Phone: 2402-2188 Fax: 2411-2323

