

AirSep® FreeStyle® 5

PORTABLE OXYGEN CONCENTRATOR

FREEDOM HAS A NEW LOOK

CAIRE introduces important updates to its high capacity portable oxygen concentrator – the AirSep® FreeStyle® 5. The new design adds features and benefits that will not only make the device more appealing to patients, but even more reliable for providers.

- **New Power Supply Connection**

Now featuring multi-directional installation with no pins to align and no release button to press, it is easier than ever to re-charge on the go.

Robust metal connection points ensure the highest durability of both the power supply and FreeStyle 5 for patients and provider alike.

- **New User-Replaceable Battery**

Swap battery packs on the go instead of waiting to re-charge an internal battery or having to add heavy external batteries. The battery pack will re-charge anytime AC or DC power is connected and is easily accessible from the front of the device.

- **Longer Battery Life**

The new battery pack provides an additional 30 minutes of duration at setting 2, providing up to 2 ½ hours with a single battery pack at setting 2.

- **New Color Scheme**

New color scheme is more appealing and non-medical in appearance. The matte finish will help providers by resisting everyday scratches.




CAIRE®
A Chart Industries Company

770.721.7759 | Toll Free 1.800.482.2473

www.CAIREMedical.com



** SEE PRODUCT WARRANTY STATEMENT FOR COMPLETE INFORMATION.

PLEASE CONSULT THE APPLICABLE PRODUCT INSTRUCTIONS FOR USE FOR PRODUCT INDICATIONS, CONTRAINDICATIONS, WARNINGS, PRECAUTIONS, AND DETAILED SAFETY INFORMATION.

© 2016 CAIRE Inc. All Rights Reserved ML-CONC0079 Rev B

PROVIDER FEATURES

- Settings 1-5 covers more patients as their disease progresses.
- Non-delivery modality eliminates expensive cylinder deliveries.
- Time-proven and trusted quality
- Standard 3-year warranty includes sieve beds

PATIENT FEATURES

- User-replaceable battery pack with up to 2½ hour battery life at setting 2
- Easily wearable with a handle, shoulder strap, or backpack deliveries.
- FAA Approved
- Fully functional on AC or DC power.
- Simple, One-Touch Operation

CLINICAL FEATURES

- UltraSense™ Technology offers sensitive breath detection
- Quickly detects and pulses to promote oxygen delivery in the earliest part of the breath
- Adjusts bolus size based on breath rate
- Up to 1000 mL/min oxygen output



SPECIFICATIONS

| | |
|---------------------------|---|
| Flow Rates | Pulse 1-5 |
| Dimensions | 10.7 x 6.6 x 4.4 in (27.2 x 16.8 x 11.2 cm) |
| Weight (w/Battery) | 6.7 lbs (3.0 kg) |
| AC Power | 100-240 VAC, 50-60Hz (2 amp max) |
| DC Power | 12-24 VDC (15 amp max) |
| Operating Temp. | 41°F to 104°F (5°C to 40°C) |
| Sound Level | 48 dB(A) @ Setting 3 |
| Operating Altitude | Up to 12,000 ft (3657.6 m) |
| O2 Concentration | 90% (+5.5% / -3%) |
| Battery Duration | Setting 1 - 4.25 hours (4:15) Setting 2 - 2.5 hours (2:30) Setting 3 - 2.0 hours (2:00) Setting 4 - 1.75 hours (1:45) Setting 5 - 1.25 hours (1:15) |
| Battery Re-Charge | 4 - 5 hours |
| Warm-Up Time | < 2 minutes |
| O2 Output | Setting 1 - 227.5 mL/min Setting 2 - 420 mL/min Setting 3 - 630 mL/min Setting 4 - 840 mL/min Setting 5 - 1000 mL/min |
| Limited Warranty | 3 years unit ** 1 year battery & power supplies |

ORDERING INFORMATION

FreeStyle 5 Includes: Portable Concentrator, Battery Pack, AC Power Supply, AC Power Cord, DC Power Supply, Carrying Case, Shoulder Strap, User Manual, Carry-all Accessory Bag

| | |
|-------------------|--------------------|
| AS077-101 | USA |
| AS077-102 | European |
| AS077-B101 | USA - 2 batteries |
| AS077-B102 | Euro - 2 batteries |
| AS077-103 | Australia |
| AS077-104 | UK |
| AS077-B104 | UK - 2 batteries |
| AS077-105 | Brazil |
| AS077-106 | Swiss |

other countries' power supplies may be available - contact customer service

ACCESSORIES

| | |
|----------------|----------------------|
| BT033-1 | Spare Battery Pack |
| MI284-1 | Backpack Harness |
| MI420-1 | Carrying Case |
| MI372-1 | Carry-All Bag |
| MI426-1 | Shoulder Strap |
| FI194-1 | Intake Filter |
| PW034-1 | DC Power Supply |
| PW035-1 | AC Power Supply (US) |
| PW035-2 | AC Power Supply (EU) |
| PW035-3 | AC Power Supply (UK) |

Cords Included

| | |
|----------------|----------------------|
| BT032-1 | Desktop Charger (US) |
| BT032-2 | Desktop Charger (EU) |
| BT032-4 | Desktop Charger (UK) |

Power Supply Included

CAIRE®

A Chart Industries Company

770.721.7759 | Toll Free 1.800.482.2473

www.CAIREMedical.com



** SEE PRODUCT WARRANTY STATEMENT FOR COMPLETE INFORMATION.

PLEASE CONSULT THE APPLICABLE PRODUCT INSTRUCTIONS FOR USE FOR PRODUCT INDICATIONS, CONTRAINDICATIONS, WARNINGS, PRECAUTIONS, AND DETAILED SAFETY INFORMATION.

© 2016 CAIRE Inc. All Rights Reserved ML-CONC0079 Rev B



Authorized Distributor:

The HomeCare Medical Ltd.
Units 1-3, 4/F, Wing Ming Industrial Centre,
15 Cheung Yue Street, Lai Chi Kok,
Kowloon, Hong Kong
Tel: 2402-2188 Fax: 2411-2323
Website: www.homecare-medical.com