

sleepcube

STANDARD

...feel refreshed!





The SleepCube Standard is the first in the new SleepCube range of Sleep therapy devices from DeVilbiss Healthcare.

The SleepCube standard is a fixed pressure CPAP designed for the effective treatment of OSA. It is equally suitable for use in the home or the hospital and its robust ergonomic design and small size makes it the ideal travelling companion. The SleepCube Standards' multi-voltage

capability means that it can be used virtually anywhere in the world and can also be operated from a 12 V DC power source such as those found in a motor home, recreational vehicle or boat.

The optional heated humidifier integrates seamlessly into the system and the SleepCube standard is packed full of all the features needed for a restful and refreshing nights sleep.



SPECIFICATIONS

Dimensions (H x W x L)

10.8 x 16.5 x 17.8 cm

Weight

1.18 kg

CPAP Pressure Range

3 - 20 cm H₂O

Operating

5° C to 40° C; 0 % to 95 % R.H.
non-condensing

Atmospheric Conditions

Sea level to 9,000 feet

Storage & Transportation

- 40° C to + 70° C; 0 % to 95 % R.H.
non-condensing

Sound Level

< 26 dB(A)

Maximum Power Consumption

65 W from AC power source
(PAP device only)

Electrical Requirements

100 - 240 VAC, 50/60 Hz



The simple LED push button display, makes the SleepCube easy to use, even in the dark



Back of the SleepCube – showing filter, patient tubing outlet and download outlet



SleepCube with integrated heated humidifier



The heated humidifier is easy to take apart and clean



DeVilbiss Healthcare Ltd.
Unit 3, Bloomfield Park, Bloomfield Road
Tipton, West Midlands DY4 9AP
United Kingdom
Phone
Fax + 44 (0) 121 521 3141
www.devilbisshc.com
enquiries@devilbisshc.com



Authorized Distributor:
The HomeCare Medical Ltd.
Units 1-3, 4/F, Wing Ming Industrial Centre,
15 Cheung Yue Street, Lai Chi Kok,
Kowloon, Hong Kong
Tel: 2402-2188 Fax: 2411-2323
Website: www.homecare-medical.com