

## Pulmo-Aide<sup>®</sup> Compressor Nebulizer System

### The Gold Standard for High Performance Medication Delivery



Pulmo-Aide Compressor  
Nebulizer System  
Model 5650D

The DeVilbiss Pulmo-Aide has been providing **effective aerosol treatments for more than thirty years**. Modern styling, durability and ease of use make this unit **perfect for in-home asthma treatments**. The Pulmo-Aide comes ready to use with all the attachments - filter, tubing, and nebulizer.

The **Pulmo-Aide is the only approved compressor for use with TOBI** (tobramycin solution for inhalation). The DeVilbiss Pulmo-Aide is unmatched in quality, value and performance.

### Performance Features

- Maximum compressor pressure – 35 psig or greater
- Compressor free air flow – 9 lpm or greater
- 1-year limited warranty
- Longest lasting compressor on the market with an average product life of 10 years
- Unmatched quality and durability - Great for nursing homes and long-term users

### Design Features

- Metal frame motor construction
- Unique 4-pole motor design extends life of the unit
- Large built-in compartment for storing nebulizers and tubing
- Built-in nebulizer holder is a convenient place to store your nebulizer if your treatment is interrupted.

## Parts and Accessories

### DEVILBISS PART #

### PRODUCT DESCRIPTION

800D	Reusable Nebulizer – blue (1 each)
800C	Reusable Nebulizer – purple/green (1 each)
4650D-621	Disposable Nebulizer (50 pack)
AAM12	Adult Aerosol Mask (12 pack)
AAM20	Adult Aerosol Mask (50 pack)
PAM10	Pediatric Aerosol Mask (50 pack)
DLI050	Pediatric Children's Character Mask (50 pack)
5650D-601	Compressor Filter (package of 25)
5650D-602	Compressor Filter (package of 5)
5650D-612	Filter Kits (6 each – caps and filters)
5650D-611	Gray Hose Tubing Connector (package of 6)

## 800 Reusable Jet Nebulizer

- MMAD 3.6 – 4 microns
- Approximately 75% of particles are below 5 microns in size
- Reduces treatment time nearly in half – 3cc treatment in 7 – 9 minutes
- 6cc capacity
- Dishwasher safe to promote frequent cleaning
- Ergonomic design fits comfortably in the hand – even with limited functionality
- 4 components for ease of assembly/disassembly
- Comfortable mouthpiece design
- Standard 22mm fitting, no adapters required
- Can be used inline with a ventilator
- 1 year warranty
- Tested for use with Albuterol, Pulmicort (budesonidel, Ipratropium bromide)
- Suitable for use with most common liquid medications

### Technical Specifications

<b>Built to meet the requirements of HCPCS code:</b> E0570		<b>Filter Replacement:</b>	No tools required
<b>Product Dimensions:</b>	10.5H X 6.3W X 10.1D	<b>Warranty:</b>	1 Year
<b>Product Weight:</b>	7.1 lbs.	<b>Safety Approvals:</b>	UL1431
<b>Maximum Air Pressure:</b>	30+ psig	<b>Motor/Pump Assembly:</b>	Diaphragm, Primarily Metal Parts
<b>Maximum Liter Flow:</b>	9+ lpm	<b>Neb Storage:</b>	Yes
<b>MMAD:</b>	5 microns or less	<b>Master Carton Quantity:</b>	4
<b>Sound:</b>	51 dBA		



DeVilbiss Healthcare | 100 DeVilbiss Drive | Somerset, PA 15501 USA  
800.338.1988 | 814.443.4881 | [www.DeVilbissHealthcare.com](http://www.DeVilbissHealthcare.com)



Authorized Distributor:  
The HomeCare Medical Ltd.

Units 1-3, 4/F, Wing Ming Industrial Centre,  
15 Cheung Yue Street, Lai Chi Kok,  
Kowloon, Hong Kong

Tel: 2402-2188

Fax: 2411-2323

Website: [www.homecare-medical.com](http://www.homecare-medical.com)