

HbO-Smart




The new HbO-Smart device is the latest finger tip pulse oximeter from Drive DeVilbiss International intended for use in homes or hospitals for non-invasive measurement of oxygen saturation, pulse rate and perfusion index.

Including adjustable settings for user alerts and four optional screen displays this new stylish device has everything you would expect from a modern day spot checking oximeter which can be used for both children and adults.

TECHNICAL SPECIFICATIONS

Display Type	OLED
Dimensions	65mm x 35mm x 30mm
Oxygen Saturation	
Range	70 - 99%
Accuracy	80 - 99% / ±2%
Resolution	1%
Pulse Rate	
Range	30 bpm - 240 bpm
Accuracy	± 2 %
Resolution	1 bpm
Perfusion Index	0.3 - 20%
Battery	2 x AAA
Battery Life	25 hours
Power Consumption	< 30mA



SpO₂, Pulse Rate and Pulse Rate Bar, Plethysmograph Waveform, Perfusion Index and Battery Status Indicator

The HbO-Smart is ready to use as it includes a lanyard, a 5 language user manual and 2 x AAA batteries.



Drive DeVilbiss International

DeVilbiss Healthcare GmbH
Kamenzer Straße 3
68309 Mannheim, Germany

UK +44 (0) 121 / 521 31 40
DE +49 (0) 621 / 1 78 98 0
kontakt@devilbisshc.com

www.drivedevilbiss-int.com

