

Founded in 1993, HomeCare Medical has been offering highly efficient products and services to patients and medical units. With close connection to hospitals and clinics, we always ensure accurate, effective and safe treatments to users in accordance with the doctor's prescription.

HomeCare Medical has been distributing the DeVilbiss Healthcare (USA) products exclusively in Hong Kong for over 20 years. The U.S. brand is one of the most commonly prescribed brands worldwide when treating various kinds of respiratory disease, with products being FDA-approved.

DeVilbiss Healthcare  
Certificate of Sole-distributor  
Authorization

Sleep Therapy on ALL CPAP  
Oxygen Therapy  
Aerosol Therapy  
Suction Therapy



DeVilbiss Healthcare awards to HomeCare Medical



Lai Chi Kok Headquarter

Units 1-3, 4/F, Wing Ming Industrial Centre, 15 Cheung Yue Street, Lai Chi Kok  
Tel : (852) 2402-2188



Kowloon Tsim Sha Tsui Showroom

Rm B, 10/F, Season Commercial Building, 3-3A Humphreys Avenue, Tsim Sha Tsui  
Tel : (852) 2992-0226



Hong Kong Island Causeway Bay Showroom

Rm 1206, 12/F, Capitol Centre, 5-19 Jardine's Bazaar, Causeway Bay  
Tel : (852) 2882-7922

www.homecare-medical.com

# SUCTION & NEBULIZER



Hong Kong Authorized Sole Distributor

www.homecare-medical.com

(852) 2402-2188

## VacuAide QSU

### Portable Suction Unit

7314P-U



With a simple on/off operation and large, rotating pressure adjustment gauge, this ultra-quiet unit allows the required suction pressure to be easily set and viewed.

- Features low noise level without compromising suction performance
- Accessories store and transport easily in accompanying carry bag
- ISO 10079-1 conformity

The AirForce Mini is a battery-operated portable inhalation therapy system. Whether in the office, at school or on the road, the AirForce Mini gives the user complete independence.

Designed using modern day compressor-technology, the AirForce Mini produces a high percentage of respirable aerosol particles suitable for treating respiratory conditions of the upper and lower airways. Exceptional performance and value; the AirForce Mini is the ideal compact solution for aerosol therapy which is suitable for both adults and children.

- One of the smallest nebulizer systems in the market
- Low noise 45dBA
- Integrated battery with LED-status display
- Stylish design
- User friendly and easy to clean
- Lightweight

Dimensions :	20.3 cm (8") x 22 cm (9") x 21.1 cm (8")
Vacuum Range :	50 - 550 mmHg $\pm$ 10%
Air Flow @ Pump Inlet :	27 l/min
Disposable Collection Container Capacity :	800ml (cc)
Reusable Collection Container Capacity :	1200ml (cc) also available
Mode of Operation :	intermittent 30 mins. on; 30 mins. off
Power Supply :	100 - 240V, 50/60 Hz
Total Weight (incl. container) :	3.0 kg (6.6 lb)

## AirForce Mini

### Compact Compressor Nebulizer

MQ6010



## AirForce One

### Compressor Nebulizer

MQ5900-UK



The AirForce One compact compressor nebulizer offers a cost-effective and reliable form of aerosol therapy. The easy-to-operate and durable unit can be used on adult and pediatric patients

- Convenient, built-in carry handle
- Comes with necessary accessories and carry bag for immediate use

Dimensions :	11.1cm (4.4") x 6.8cm (2.7") x 4.5cm (1.8")
Particle Size (MMAD) :	4.12 $\mu$ m
Air Flow :	0.15 l/min
Nebulization Rate :	0.25 ml/min
Power Supply :	230V, 50Hz
Total Weight :	0.217 kg (0.5 lb)

Dimensions :	14 cm (6") x 19 cm (7") x 9 cm (4")
Particle Size (MMAD) :	3.16 $\mu$ m
Noise Level :	52dBA
Air Flow :	5.31 l/min
Nebulization Rate :	0.4 ml/min
Power Supply :	230V, 50Hz
Total Weight :	1.5 kg (3.3 lb)

## PulmoAide

### Compressor Nebulizer

5650I



Modern styling, durability and ease of use make this unit perfect for in-home asthma treatments.

- Portable, lightweight and robust design
- Comes with disposable nebulizer
- Built-in storage compartment to hold nebulizer and tubing until needed again
- Integrated handle for easy carrying
- Easy set-up within seconds the PulmoAide is ready
- 75 % of the medication is deposited in the lungs with an MMAD of 3.6 - 4.0 mm
- High Aerosol output to reduce treatment time

Dimensions :	17 cm (6.5") x 26 cm (10.1") x 27 cm (10.5")
Nebulization Rate :	>0.15 ml/min
Compressor Free Air Flow :	12 l/min
Supplied Nebulizer MMAD :	<5 $\mu$ m
Power Supply :	230V, 50Hz
Total Weight :	3.2 kg (7.1 lb)

## UltraNeb

### Ultrasonic Nebulizer

U3000



This high-output, ultrasonic nebulizer is optimal for humidification, sputum removal and aerosol therapy.

- Air flow output and nebulization rates have four adjustable settings and timer function from 0-60 minutes
- Integrated alarm functions alert when chamber is not attached or empty, medication time is complete, or device fails
- Suitable for most sterile water systems
- Compatible with various conveyor systems

Dimensions :	20.5 cm (8") x 32 cm (12") x 19 cm (7")
Nebulization Rate :	3 ml/min
Air Flow :	20 l/min max
Particle Size (MMAD) :	<4 $\mu$ m
Ultrasonic Frequency :	1.68 MHz $\pm$ 5%
Power Supply :	230V, 50Hz
Total Weight :	3.5 kg (7.7 lb)